

Certified Blue TOC Vials

FACT SHEET

Certified <10 ppb TOC and DUCT vials with clean packaging for controlled environments

Veolia—the global leader in low-level TOC standards and TOC and conductivity measurement technology through its Sievers* product line—offers Certified Blue packaging for its suite of pre-cleaned TOC and Dual Use Conductivity & TOC (DUCT) vials.

Veolia understands the need to control the cleanliness of materials in manufacturing and laboratory spaces, including paper and cardboard products that are potential sources of microbial contamination.

With Sievers Certified Blue vial products, conventional corrugated cardboard boxes are replaced with recyclable, blue polypropylene boxes. Paper certificates of conformance and paper labels are replaced with plastic labels and barcodes for easy access to the electronic certificates. The Clean Package design, triple-bagged TOC 12-pack or DUCT 30-pack, is ideal for cleanrooms.

Sievers Certified Blue TOC vials are designed for use in critical pharmaceutical applications such as USP, EP, JP, and KP Water testing. Our certification to <10 ppb



Available in 9 configurations

TOC (and no ionic leaching for DUCT vials) or <35 ppb TOC for pre-acidified vials gives you confidence in the quality of TOC and conductivity measurements. Additionally, the vials are scrupulously cleaned and assembled in an ISO Class 7 environment to minimize particulates. Sievers Certified Blue TOC vials are compatible with most TOC analyzer brands, including Sievers, Shimadzu, OI Analytical, and Teledyne-Tekmar.

	Product Name	Part Number	Count	Product Details	Packaging
DUCT Vials 30-Pack	Certified Blue DUCT Vials	HMI 77510-01	30 vials	DUCT vial, <10 ppb TOC, No ionic leaching	Polypropylene Box, Single Bag
	Certified Blue DUCT CV Kit	HMI 77520-01	30 vials 40 swabs	Swab: <50 ppb TOC No ionic leaching DUCT vial: <10 ppb TOC	Polypropylene Box, Single Bag, (2) 20 pack Swabs
	Certified Blue DUCT Clean Package	HMI 77530-01	30 vials	No ionic leaching DUCT vial, <10 ppb TOC	Polypropylene Box, Triple Bagged
TOC Vials 72-Pack	Certified Blue TOC Vials	HMI 90610-01	72 vials	40 mL, TOC vial <10 ppb TOC	Polypropylene Box, Single Bag
	Certified Blue Large CV Kit	HMI 90620-01	72 vials 80 swabs	40 mL TOC vial: <10 ppb TOC, Swab: <50 ppb TOC	Polypropylene Box, Single Bag, (4) 20 pack Swabs
	Certified Blue Pre-acidified Vials	HMI 90690-01	72 vials	40 mL TOC vial, 20 uLH ₃ PO ₄ < 35 ppb TOC	Polypropylene Box, Cleanroom Sorbent Pad, Single Bag
TOC Vials 12-Pack	Certified Blue 12-Pack	HMI 91210-01	12 vials	40 mL TOC vial, <10 ppb TOC	Polypropylene Box, Single Bag
	Certified Blue Small CV Kit	HMI 91220-01	12 vials 20 swabs	40 mL TOC vial: <10 ppb TOC, Swab: <50 ppb TOC	Polypropylene Box, Single Bag, (1) 20 pack Swabs
	Certified Blue Clean Package	HMI 91230-01	12 vials	40 mL TOC vial, <10 ppb TOC	Polypropylene Box, Triple Bagged

Veolia Water Technologies

Please contact us via:

www.veoliawatertechnologies.com